



# Rogozinski Lamina Spreader


Designed by Chaim Rogozinski, MD

*Self-retaining and self-leveling lamina spreader that captures the spinous processes, thereby helping to maintain interlaminar retraction*



#### PRODUCT NO'S:

4275 [Standard]
Overall Length: 8" (20,3 cm)
Leg Depth: 3" (7,6 cm)
Blade Width: .75" (1,9 cm)
Blade Height: .5" (1,3 cm)
4275-01 [Small]
Overall Length: 8" (20,3 cm)
Leg Depth: 2" (5,1 cm)
Blade Width: .75" (1,9 cm)
Blade Height: .5" (1,3 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362