


Rogozinski Soft Tissue Retractor

Designed by Chaim Rogozinski, MD

Self-leveling retractor that helps lessen tissue movement underneath the prongs, thereby helping to maximize exposure



PRODUCT NO'S:	
4276 [Standard]	
Overall Length: 8" (20,3 cm)	
Leg Depth: 2.5" (6,4 cm)	
Blade Width: .75" (1,9 cm)	
Blade Height: 1.5" (3,8 cm)	
4276-01 [Small]	
Overall Length: 8" (20,3 cm)	
Leg Depth: 2" (5,1 cm)	
Blade Width: .75" (1,9 cm)	
Blade Height: 1.5" (3,8 cm)	



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362