

Sarraf TiN Coated Cement Forceps

Designed by Khaled M. Sarraf, MD

TiN coating helps reduce metal transfer and protects the prosthetic implant surface when in contact

Aids with the removal of excess extruded cement during joint arthroplasty.



Ultra hard titanium nitride coating helps to extend forceps life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion, while helping to eliminate metal transfer and protect the implant surface.

PRODUCT NO'S:

5039 [Straight]

Overall Length: 6" (15,2 cm)

5041 [Angled]

Overall Length: 6.125" (15,6 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362