

Sarraf TiN Coated Cement Forceps

Designed by Khaled M. Sarraf, MD

TiN coating helps reduce metal transfer and protects the prosthetic implant surface when in contact

Aids with the removal of excess extruded cement during joint arthroplasty, and can be used in conjunction with any Innomed TiN coated instrument.



Ultra hard titanium nitride coating helps to extend forceps life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion, while helping to eliminate metal transfer and protect the implant surface.

PRODUCT NO'S:

5039 [Straight]
Overall Length: 6" (15,2 cm)

5041 [Angled]
Overall Length: 6.125" (15,6 cm)

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362