

Shoulder Surgery Retractor System

Developed in collaboration with Mayo Clinic.

PRODUCT NO'S:

1251-00 [Complete System]

Included in Set/Available Individually:

1252-N [Modified Thin Glenoid Retractor–Narrow]

Two included in set; one with this product number
Overall Length: 11.875" (30,2 cm)
Blade Width: 15 mm

1252-W [Modified Thin Glenoid Retractor–Wide]

Two included in set; one with this product number
Overall Length: 11.875" (30,2 cm)
Blade Width: 23 mm

1253 [Right Angle Hohmann Retractor]

Overall Length: 8.125" (20,6 cm)
Depth from Bend: 4.25" (10,8 cm)
Blade Width: 16 mm

1254 [Modified Fukuda Retractor]

Overall Length: 8.625" (21,9 cm)
Depth: 2.75" (7 cm)
Blade Width: 39 mm

1255-L [Brown Deltoid/Richardson Retractor–Large]

Overall Length: 10.5" (26,7 cm)
Depth: 2.5" (6,4 cm)
Blade Width: 60 mm

1255-S [Brown Deltoid/Richardson Retractor–Small]

Overall Length: 10.5" (26,7 cm)
Depth: 2.5" (6,4 cm)
Blade Width: 44 mm

1256 [Modified Darrach Retractor, Straight–Narrow]

Overall Length: 10.25" (26 cm)
Blade Width: 12,7 mm

1257 [Modified Darrach Retractor, Straight–Wide]

Overall Length: 10.25" (26 cm)
Blade Width: 19 mm

1258 [Modified Darrach Retractor, Bent–Narrow]

Overall Length: 10.75" (27,3 cm)
Blade Width: 12,7 mm

1259 [Modified Darrach Retractor, Bent–Wide]

Overall Length: 10.75" (27,3 cm)
Blade Width: 19 mm

1260 [Soft Tissue Shoulder Retractor]

Overall Length: 10" (25,4 cm)
Depth from Bend: 3" (7,6 cm)
Blade Width: 19 mm

1261 [Glenoid Access Retractor]

Overall Length: 13.5" (34,3 cm)
Access Hole Internal Diameter: 36 mm X 30 mm
Depth of Prongs: 8.5 mm



System includes two of each size of the Modified Thin Glenoid Retractors, and one of each of the other retractors.

Modified Thin Glenoid Retractor–Narrow & Wide



Right Angle Hohmann Retractor



Modified Fukuda Retractor



Brown Deltoid/Richardson Retractor–Large & Small



Modified Darrach Retractor, Straight–Narrow & Wide



Modified Darrach Retractor, Bent–Narrow & Wide



Soft Tissue Shoulder Retractor



Glenoid Access Retractor



New!

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2020 Innomed, Inc.