


Slavitt Phalangeal Forceps

Designed by Jerome Slavitt, DPM

Designed to enable the surgeon to provide joint distraction and stability during prosthesis placement at the base of the proximal phalanx of the lesser digits



Helps to distract the joint and hold the bone, allowing easier access to the base. Can also be used for digital fusions to hold bones better for drilling and cutting applications.

PRODUCT NO:
1163
Overall Length: 6" (15,2 cm)
Clamp Internal Opening Diameter: 4 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.