

Anderson Talar Neck Osteotome

Designed by John Anderson, MD


Designed to help improve range of motion and reduce pain caused by anterior bony impingement of the ankle by removing osteophytes from the anterior talar neck and the anterior distal tibia



Three Widths Available

PRODUCT NO'S:

5075 Osteotome Width: 17 mm Overall Length: 9.875" (25,1 cm) Handle Length: 4.5" (11,4 cm)	5075-50 Osteotome Width: 12.7 mm Overall Length: 9.875" (25,1 cm) Handle Length: 4.5" (11,4 cm)	5075-75 Osteotome Width: 9.5 mm Overall Length: 9.875" (25,1 cm) Handle Length: 4.5" (11,4 cm)
---	--	---



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.