


# Unger Anterior Total Hip Instruments

Designed by Anthony Unger, MD

Universal instrument system specifically designed for Direct Anterior approach THR

## PRODUCT NO.'S:

1	3001 [Unger Wide Hohmann-Single Prong] Blade Width: 45 mm Blade Depth: 5.45" (13,8 cm) Overall Length: 13.5" (34,3 cm)
2	3008 [Unger Wide Hohmann-Double Prong] Blade Width: 45 mm Blade Depth: 5.45" (13,8 cm) Overall Length: 13.5" (34,3 cm)
3	3002 [Unger Narrow Hohmann] Blade Width: 34 mm Blade Depth: 4" (10,2 cm) Overall Length: 13" (33 cm)
4	3003 [Unger Blunt Narrow Cobra] Blade Width: 34 mm Blade Width at Tip: 12 mm Blade Depth: 5.25" (13,3 cm) Overall Length: 14.5" (36,9 cm)
5	3004 [Unger Canal Finder Rasp-Straight] Overall Length: 11" (27,9 cm) Handle Length: 5" (12,7 cm)
6	3004-01 [Unger Canal Finder Rasp-Curved] Overall Length: 11" (27,9 cm) Handle Length: 5" (12,7 cm)
6s	3004-02 [Unger Canal Finder Rasp-Curved with Smooth Proximal] Overall Length: 11" (27,9 cm) Handle Length: 5" (12,7 cm)
7	3005-R [Unger Box Osteotome-Right] Overall Length: 12" (30,5 cm)
8	3005-L [Unger Box Osteotome-Left] Overall Length: 12" (30,5 cm)
9	3006 [Unger Femoral Neck Elevator] Blade Width at Widest: 25 mm Overall Length: 13" (33 cm) Handle Length: 9" (22,9 cm)
10	3007 [Unger Soft Tissue Protector] Blade Width: 50 mm Blade Depth: 1.75" (4,4 cm) Overall Length: 10.125" (25,7 cm)
11	3009-R [Unger Offset Narrow Hohmann-Right] Blade Width: 34 mm Blade Depth: 4" (10,2 cm) Overall Length: 13" (33 cm)
12	3009-L [Unger Offset Narrow Hohmann-Left] Blade Width: 34 mm Blade Depth: 4" (10,2 cm) Overall Length: 13" (33 cm)
<b>Optional Instruments:</b>	
13	3006-01 [Femoral Neck Elevator-Long Prong] Blade Width at Widest: 25 mm Overall Length: 13.4" (34 cm) Handle Length: 9" (22,9 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.