

# Universal Screwdriver Set

*Helps eliminate the opening of multiple sterile packs when a specific size or style of screwdriver is needed*

Helpful during revision total joint surgery where screws have been used, removal of bone plates, fracture fixation screws or bone graft screws. Set consists of: 7 (seven) double ended screwdriver bits – small & large single slot, cross & cruciate, 3.5 mm & 4.5 mm hex, small & large phillips, small, medium, & large star – a handle which accommodates any of the above bits, and a sterilization case.



Set in Storage Case



PRODUCT NO'S:
5195 [Complete Set with Case] <i>Also sold individually</i>
5195-01 [Handle]
5195-02 [Straight (single slot)] <i>Large: 7 x 1.5 mm, Small: 5 x 1 mm</i>
5195-03 [Cross/Cruciate] <i>Large: 7 mm, Small: 6 mm</i>
5195-04 [Hex] <i>Large: 4.5 mm, Small: 3.5 mm</i>
5195-05 [Phillips] <i>Large: 4 mm, Small: 3.5 mm</i>
5195-08 [Small Star: #6 & #8]
5195-06 [Medium Star: #10 & #15]
5195-07 [Large Star: #20 & #25]



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**