



Wagner Osteotome Handle

Handle is designed for easier gripping, rotational control, and use with a mallet with a standard 1/4" Lambotte osteotome



PRODUCT NO'S:
5348 [Handle Only] Overall Length: 5.5" (14 cm)
5348-01 [1/4" Osteotome Only] Overall Length: 8.875" (22,5 cm)

Handle designed by
Russell Wagner, MD



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362