

Weatherly Mini-Deltoid Retractors

Designed by Wallace Weatherly, MD

Designed for the retraction of the deltoid in a mini-open mid-deltoid splitting approach to rotator cuff surgery, the offset handle helps allow clear visualization of the surgical field, and the ergonomic non-slip handle surface helps prevent fatigue in the operative team



PRODUCT NO'S:
5110-L [Large] Overall Length: 12.75" (32,4 cm) Depth from Bend: 4.5" (11,4 cm) Blade Dimensions: 40 mm x 90 mm
5110-M [Medium] Overall Length: 11" (27,9 cm) Depth from Bend: 3" (7,6 cm) Blade Dimensions: 40 mm x 55 mm
5110-S [Small] Overall Length: 10.5" (26,7 cm) Depth from Bend: 2.5" (6,4 cm) Blade Dimensions: 40 mm x 40 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362