

Wiater Shoulder Bone Hook


Designed by J. Michael Wiater, MD, FAOS, FAOA

Large bone hook designed to retract the proximal humerus posteriorly to help provide exposure for glenoid reaming during open shoulder procedures

Also useful for other large joint procedures.



PRODUCT NO:
5079
Curve Diameter: 2.68" (68 mm)
Overall Length: 10" (25,4 cm)
Handle Width: 4" (10,2 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED