


Wilson Trigger Finger Retractor

Designed by Ralph V. Wilson, MD

Helps provide improved exposure, which frees up the assistant and gives easier and better visualization of the tendon sheath in a trigger finger procedure

Also useful in other small incision surgeries. Can use with Morton's neuroma and hand arthroplasty procedures.



PRODUCT NO:

1884

Overall Length: 4.25" (10.8 cm)

Blades: 6.5 mm Wide x 10 mm Deep

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.