



Extra Large Hip Retractors

Designed by Wayne M. Goldstein, MD



For hip surgery with large patients, and when extra large instruments are desired for increased leverage and depth


Extra Leverage Femoral Neck Elevators

PRODUCT NO'S:	
7650 [Standard]	Overall Length: 18.25" (46,4 cm) Handle Length: 9.25" (23,5 cm) Blade Width: 38 mm
7650-02 [Short Handle]	Overall Length: 15.25" (38,8 cm) Handle Length: 6.25" (15,9 cm) Blade Width: 38 mm



Infero-posterior Acetabular Capsule Retractors

PRODUCT NO'S:	
7620-01 [Right]	Overall Length: 12" (30,5 cm) Handle-to-Bend Length: 6" (15,2 cm)
7620-02 [Left]	Overall Length: 12" (30,5 cm) Handle-to-Bend Length: 6" (15,2 cm)



Extra Leverage Proximal Femoral Elevator

PRODUCT NO:	
7640	Overall Length: 17.5" (44,5 cm) Handle Length: 13" (33 cm) Blade Width at Widest: 63 mm



Large Cobra Retractors

PRODUCT NO'S:	
7630-01 [Standard]	Overall Length: 17.5" (44,5 cm) Handle Length: 14" (35,6 cm)
7630-02 [Wide]	Overall Length: 17.5" (44,5 cm) Handle Length: 14" (35,6 cm)



Modular Weights

Weights can be used to help hold the retractors in place.

PRODUCT NO'S:	
3430-01	1.5 lbs. (.68 kg)
3430-02	2.0 lbs. (.91 kg)
3430-03	2.5 lbs. (1.13 kg) with attaching hook



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362